



BPB MASK



FEATURES:

- Single use surgical mask type IIR
- MADE IN ITALY
- BPB MASK guarantees nose and mouth basic protection from external agents such as dust, pollen and bacteria.
- Made with layers of SMS (Thermobonded Nonwoven 100% PP Hydrophobic), equipped with adjustable nose clip and ear elastics
- Bacterial Filtration Efficiency BFE 99%
- **MEDICAL DEVICE - EC conformity according to MDD 93/42/EC - Class I non sterile**
- Sizes: 17,5 x 9,5 cm



PRODUCT CODE

BPB MASK-001



MADE IN ITALY



BPB MASK CATEN
Rev.02 15/06/2020

Manufacturer:

BIOPSYBELL S.R.L.

Via Aldo Manuzio 24 - 41037 Mirandola (MO) Italy
Phone +39 0535 27850 – Fax +39 0535 33526
C.F./P.Iva 02615000367 - www.biopsybell.it

Contact us for more information:

international1@biopsybell.it

italia@biopsybell.it